

Medtronic's CD Horizon® Legacy™ Peek Rod System

Last Updated Tuesday, 15 May 2007

See an innovative lumbar fusion using Medtronic's CD HORIZON® LEGACY™ PEEK Rod System. Join Todd Lanman, M.D. and Russ Nockels, M.D. as they perform this procedure during a live surgical broadcast at 6:00 PM CDT (23:00 UTC) on OR-Live.com

The CD HORIZON LEGACY PEEK Rod System uses rods made of polyetheretherketone, a semi-crystalline thermoplastic polymer that meets all biocompatibility requirements and has a long history of use in surgical implants. When combined with CD HORIZON® Spinal System-based pedicle screws and minimally invasive approach options such as the MAST QUADRANT™ Retractor Set, the CD HORIZON LEGACY PEEK Rod dynamically stabilizes the affected segment while allowing for a streamlined surgical technique that preserves much of the natural anatomy.